

Devyser Complete

CE IVD	Yes
Validated PCR instruments	ABI/PE TC 9600 ABI/PE TC 9700 ABI/PE TC 2720 Eppendorf Mastercycler MJ Research (BioRad) PTC-200 series
Validated genetic analyzers	ABI Prism® 310/3100/3130/3700 MegaBACE
Validated polymers	POP-4 / POP-6 / POP-7
Validated sample preps	QIAamp Blood DNA Mini Kit QIAamp DNA Mini Kit
Tested sample preps	Qiagen EZ1 MagNA Pure InstaGene Matrix
Reaction volume (µL)	25 µL each mix / clinical sample
No. of markers chromosomes 13, 18, 21	21 (7 each chromosome)
No. of sex chromosome markers	7 (2 sex determining in Mix 2)
Unique marker for Turner Syndrome	Yes
No. of markers total	28
Detection range	1 - 30 ng genomic DNA / PCR
Ready-To-Use for PCR	Yes
Markers designed to allow data analysis with poor quality DNA	Yes
Positive and Non-template control included	Yes (54 test kit)
Available sizes	9/54/100 test kits

Devyser AB  
Instrumentvägen 19  
SE-12653, Hägersten  
Sweden

Phone: +46-8-56215850  
Fax: +46-8-32 64 88  
www.devyser.com